

Footrest and Legrest Accessories

| Item Number | Note | Part Number | Description | Quantity Per | |
|-------------|-------|--------------------|---|--------------|----------|
| | | | | Package | Assembly |
| | G,I | 1026987 | Guard, Impact (8" L) | 1 | 1 |
| | G,H | 1143050 | Non-Removable Device, Footrest/Legrest (1820) | 1 | 1 |
| | | | NLA - Non-Removable Device, Footrest-Legrest | 1 | 1 |
| | | | Carton | 1 | 1 |
| | | 1032002 | Strap, "H" Heel | 1 | 1 |
| | | 1035902 | Strap, Calf w/ Hook and Loop (31" L) | 1 | 1 |
| | A,B,C | 8881006725-NL A | NLA - Panel, Legrest | 1 | 1 |
| | A,B,D | 8881006726 | Panel, Legrest | 1 | 1 |
| | A,B,E | 8881006727-NL A | nla - Panel, Legrest | 1 | 1 |
| | A,C,F | 8881035816-NL A | NLA - Panel, Legrest | 1 | 1 |
| | A,D,F | 8881035817-NL A | NLA - Panel, Legrest | 1 | 1 |
| | A,E,F | 8881035818-NL A | NLA - Panel, Legrest | 1 | 1 |
| | G | See Chart | Heel Loops and Attaching Hardware | 1 | 1 |
| | G | See Chart | Heel Loop Attaching Hardware | 1 | 1 |
| | G | See Chart | Toe Loops and Attaching Hardware | 1 | 1 |
| | G | See Chart | Toe Loop Attaching Hardware | 1 | 1 |

NOTE: A - Specify color when ordering. See page titled, "Frame Finishes, Touch-Up Paint and Upholstery Options."
 B - PW93E.
 C - 16" Wide.
 D - 18" Wide.
 E - 20" Wide.
 F - P904A and PAL4A.
 G - Not Shown.
 H - Includes items Device, Non-Removable and Carton.
 I - For 70 deg. Tapered Footrest.

Heel Loops and Attaching Hardware

| Footplate Type | Heel Loop Single # | Heel Loop Pair # | Heel Loop w/ Ankle Strap Single # Rt. | Heel Loop w/ Ankle Strap Single # Lt. | Heel Loop w/ Ankle Strap Pair # | Attaching Hardware Pair # |
|---|--------------------|------------------|---------------------------------------|---------------------------------------|---------------------------------|---------------------------|
| Composite (1651) Small 16" Wide Tapered | 1028412-NLA | 1035912-NLA | 1035908-NLA | 1035909-NLA | 1035906 | 1035931 |
| Composite (1651) Medium 16" Wide (18" Wide Tapered) | 1017621 | 1035917-NLA | 1035919 | 1035920 | 1035906 | 1035931 |

Footrest and Legrest Accessories

Heel Loops and Attaching Hardware

| Footplate Type | Heel Loop Single # | Heel Loop Pair # | Heel Loop w/ Ankle Strap Single # Rt. | Heel Loop w/ Ankle Strap Single # Lt. | Heel Loop w/ Ankle Strap Pair # | Attaching Hardware Pair # |
|--|--------------------|------------------|---------------------------------------|---------------------------------------|---------------------------------|---------------------------|
| Composite (1651) Large 18-20" Wide (20" Wide Tapered) | 1017621 | 1035917-NLA | 1035932-NLA | 1035933-NLA | 1035930 | 1035931 |
| Standard Aluminum (1650) 16-20" Wide | 1017621 | 1035917-NLA | 1035919 | 1035920 | 1035918 | 1035907 |
| Deluxe Aluminum (1450) 16-20" Wide | 1017621 | 1035917-NLA | 1035919 | 1035920 | 1035918 | 1086112-NLA |
| Extra Large Aluminum (1350) 16-20" Wide | 1017621 | 1035917-NLA | 1035919 | 1035920 | 1035918 | 1055241 |

Heel Loop Attaching Hardware

| |
|---|
| Part Number |
| 1035931 |
| Kit Includes: (2) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (2) Spacer (1/4 x 11/32 x 1-3/4") (2) Hex Nut w/ Nylon Cap (1/4-20) (2) Flat Washer (1/4 x 5/8 x 1/16") |
| 1055241 |
| Kit Includes: 2 ea. Screw, Phillips Pan Head (1/4-20 x 2-1/4") 2 ea. Spacer (1/4 x 11/32 x 1-3/4") 2 ea. Locknut (1/4-20) |

Toe Loops

Footrest and Legrest Accessories

Toe Loops

| Toe Loops (1345) | Toe Loop Single # | Toe Loop Pair # | Attaching Hardware Pair # |
|---------------------|----------------------|--------------------|---------------------------------|
| All Sizes | 04572X025 | 1037087 | 1037089 |

Footrest and Legrest Accessories

Toe Loop Attaching Hardware

| |
|---|
| Part Number |
| 1037089 |
| Kit Includes: (2) Screw, Phillips Pan Head (10-32 x 7/8") (2) Washer (1/4 x 3/4 x 3/32") (2) Locknut (10-32) |