


One Piece Footboard - T1263



One Piece Footboard - T1263

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1102205	One Piece Footboard Assembly (20" W)	1	1
1	A	1102206	One Piece Footboard Assembly (22" W)	1	1
1	A	1102207	One Piece Footboard Assembly (24" W)	1	1
1	A	1102208	One Piece Footboard Assembly (26" W)	1	1
1	A	1102209	One Piece Footboard Assembly (28" W)	1	1
1	A	1102210	One Piece Footboard Assembly (30" W)	1	1
2		1005383	Screw, Socket Button Head (#10-32 x 1")	1	8
3		1108599	Plate, Lifter, Black	1	1
4		47352X004	Clamp, Black	1	4
5		1016346	Cap, End, Tube, Black (7/8")	1	2
6		1102200	Footboard (20" W)	1	1
6		1102201	Footboard (22" W)	1	1
6		1102202	Footboard (24" W)	1	1
6		1102203	Footboard (26" W)	1	1
6		1102204	Footboard (28" W)	1	1
6		1102138	Footboard (30" W)	1	1
7		1013218	Washer (#10 Flat)	1	16
8		60279X003	Locknut (#10-32)	1	8
9		1101478	Tube, Pivot/Slide, Chrome	1	2
10		1027097	Plug Button, Black (3/4")	1	2
11		47352X003	Lock, Pad Plate, Legrest, Black	1	2
12		1008432	Screw, Socket Button Head (1/4-20 x 1-1/4")	1	4
13		60275X000	Locknut (1/4-20)	1	4
14		1098801	Insert, Support, Legrest	1	2
	B	1104905	Instruction Sheet, 1-Piece Footboard	1	1

NOTE: A - Includes Items 2-14, Instruction Sheet and Cartons.
B - Not Shown.