


Crossbraces & Hardware

Removable Arm Chairs - Recliner Chairs Only

NOTE: Due to the possibility of various manufacturing dates for Tracer IV Wheelchairs built with a serial # starting with "03H" or "03I", replacement parts will need to be ordered thru an Invacare Parts Specialist at 1-800-333-6900.



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1062848-NLA	NLA - Crossbrace without Hardware, Chrome - Right (16" W)	1	1
1		1062849-NLA	NLA - Crossbrace without Hardware, Chrome - Left (16" W)	1	1
1		1062182-NLA	NLA - Crossbrace without Hardware, Chrome - Right (18" W)	1	1
1		1062183-NLA	NLA - Crossbrace without Hardware, Chrome - Left (18" W)	1	1
2		1054951-NLA	NLA - Guide, Seat	1	2
3		1053069	Insert, Seat Rail	1	2
4		1027095	Plug Button, Black (7/8")	1	2
5		1017985	Tube, Insert	1	2
6		40183X002	Screw, Phillips Pan Head (#10-32 x 1-1/8")	1	4
7		40317X001	Screw, Hex Head (5/16-24 x 2-3/4")	1	1
8		60279X003	Locknut (#10-32)	1	4
9		1052218-NLA	NLA - Seat Guide	1	2
10		1107321	Nut, Low Cap (5/16"-24)	1	1