


Permanent Arms

NOTE: Due to the possibility of various manufacturing dates for Tracer IV Wheelchairs built with a serial # starting with "03H" or "03I", replacement parts will need to be ordered thru an Invacare Parts Specialist at 1-800-333-6900.



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		40120X001	Screw, Phillips Pan Head (#8-32 x 3/8")	1	4
2	A	8881091684	Pad, Arm, Desk Length, Vinyl	1	1
2	A	8881091688	Pad, Arm, Desk Length, Lightweight Vinyl	1	1
2		1041062	Pad, Arm, Upper/Lower, Black Plastic, Desk Length	1	1
2	A	8881091682	Pad, Arm, Full Length, Vinyl	1	1
2	A	8881091686	Pad, Arm, Full Length, Lightweight Vinyl	1	1
2		1041061	Pad, Arm, Upper/Lower, Black Plastic, Full Length	1	1
3		47234X007	Guard, Clothing, Permanent Fixed Frame, Stainless Steel - Right/ Left	1	1
4		40183X002	Screw, Phillips Pan Head (#10-32 x 1-1/8")	1	2
	B	1042118-NLA	NLA - Decal, Clothing Guard - Tracer	1	1
	B	1101333-NLA	NLA - Decal, Invacare Logo (2.75)	1	1

NOTE: A - Specify Upholstery Color.
B - Not Shown.