


Forks and Hardware

NOTE: Due to the possibility of various manufacturing dates for Tracer IV Wheelchairs built with a serial # starting with "03H" or "03I", replacement parts will need to be ordered thru an Invacare Parts Specialist at 1-800-333-6900.



8" Fork and Stem

6" Fork and Stem

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1001843	Cap, Head Tube, Black	1	1
1		1131189	Package, Head Tube Cap	10	1
2		1001979	Locknut (1/2-20)	1	1
3		1021508	Washer, Nylon (5/8 x 1 x 3/32")	1	1
4		1001801	Bearing, Flange (1/2")	1	2
5	A	1037420	Fork and Stem Assembly, Black (6 x 1" Caster)	1	1
5	B	1012124	Fork and Stem, Chrome	1	1
5	C	1015789-NLA	NLA - Fork and Stem Assembly, Black	1	1
5	D	1075529	Fork and Stem, Steel, Black (8" Caster)	1	1
5	E	1013896-NLA	NLA - Fork and Stem, Chrome	1	1

NOTE: A - 6" Casters.
 B - 8 X 1" Casters - Used on Tracer, Plus, LX, IV and Recliner Before 2/7/00.
 C - 8 X 1" Casters - Used on Tracer LT and Transport.
 D - 8 X 1" Casters - Used on Tracer IV Between 2/8/00 - 9/30/03.
 E - 8 X 2" Casters.