


Footrest Support Assembly

T963A/C/S



Footrest Support Assembly T963A/C/S

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1020719	Anti-Rattle, Black	1	1
2	A	60279X003	Locknut (#10-32)	1	1
2	B	1025735	Locknut (#10-32)	1	1
3		47565X057	Bushing (7/32 x 3/8 x 11/32")	1	1
4	A	40163X008	Screw, Phillips Pan Head (#10-32 x 11/16")	1	1
4	B	1025716	Screw, Phillips Pan Head (#10-32 x 11/16")	1	1
5	A	40163X002	Screw, Phillips Round Head (#10-32 x 1/2")	1	1
5	B	1026114	Screw, Phillips Pan Head (#10-32 x 1/2")	1	1
6	A,D	1035539	Handle, Release, Black - Right	1	1
6	A,D	1035540	Handle, Release, Black - Left	1	1
6	B,D	1035539	Handle, Release, Black - Right	1	1
6	B,D	1035540	Handle, Release, Black - Left	1	1
	F	1087698-NLA	NLA - Handle, Release, Chrome - Right	1	1
	F	1087699-NLA	NLA - Handle, Release, Chrome - Left	1	1
7	A	1008935	Spring, Release - Right	1	1
7	A	1008936	Spring, Release - Left	1	1
7	B	1028707	Spring, Release - Right	1	1
7	B	1028708	Spring, Release - Left	1	1
8	A	67005X000	Rivnut (#10-32)	1	1
8	B	1028846-NLA	NLA - Rivnut (#10-32)	1	1
9	A	11232X009	Bolt, Lug (5/16-24 x 1-3/8")	1	1
9	B	1026030	Bolt, Lug (5/16-24 x 1-3/8")	1	1
10	A	1107321	Nut, Low Cap (5/16"-24)	1	1
10	B	1025730	Locknut (5/16-24)	1	1
11	A	1111153-NLA	NLA - Support, Footrest, Black - Right	1	1
11	A	1111154-NLA	NLA - Support, Footrest, Black - Left	1	1
11	B,C	8881111151-NL A	NLA - Support, Footrest, Color - Right	1	1
11	B,C	8881111152-NL A	NLA - Support, Footrest, Color - Left	1	1
12	E	See Chart	Footrest Support Assembly - T963A/C/S	1	1

NOTE: A - For use on chrome chairs.
B - For use on color chairs.
C - Specify frame color when ordering, see category titled, "Frame Finishes, Upholstery Colors, Decals and Packaging."
D - When front riggings are ordered for a 9000XT Chair with clamp-on wheel locks, the front rigging should be ordered through Invacare Specials due to the different release handle style.
E - Includes items 1 thru 11.
F - Not shown. Order this item when ordering chair with Cyclical Lever Drive only.

Footrest Support Assembly - T963A/C/S

Prior Model Number	Current Model Number	Chair Width	Right Side	Left Side
N/A	T963 A/C Color	14-22"	8881111315 See Note(s) A	8881111316 See Note(s) A
N/A	T963 A/C Black	14-22"	1111319	1111320

Footrest Support Assembly

T963A/C/S

NOTE: A - Specify Frame Color; see section titled, "Frame Finishes, Upholstery Colors, Decals, and Packaging".