

Wheels - Non-Projection Handrims and Attaching Hardware

Handrims, Chrome and Aluminum

Model Number	Wheel Type	Chrome Handrim	Aluminum Handrim	Quantity	Attaching Points	Attaching Hardware
2244BK	24" Composite Pneumatic	1054966	1015406	1	7	1025153
U2244BK	24" Composite Urethane	1054966	1015406	1	7	1025153
4244BK	24" Composite Pneumatic or Flat Free	1054966	1015406	1	7	1025153
U2222	24" Composite Urethane	1055272-NLA See Note(s) A	N/A	1	6	1042074
2222BK	24" Composite DMT	1055272-NLA See Note(s) A	N/A	1	6	1042074
1764BK	22" Composite Pneumatic	1020000	1020018	1	7	1025153
4764 BK	22" Composite Pneumatic or Flat Free	1020000	1020018	1	7	1025153
U1764BK	22" Composite Urethane	1020000	1020018	1	7	1025153
1763BK	20" Composite Urethane	1019999	1020017	1	7	1025153
4763BK	20" Composite Pneumatic or Flat Free	1019999	1020017	1	7	1025153
U1763BK	20" Composite Urethane	1019999	1020017	1	7	1025153

NOTE: A - Chrome Handrim part #1055272 replaced with Clear Anodized Handrim part #1183681.

NOTE: Kits Contain Hardware for One Handrim

Handrims, Composite

Model Number	Wheel Type	Handrim	Quantity	Attaching Points	Attaching Hardware
2244BK, 4244BK, U2244BK	24" Composite Pneumatic or Flat Free	1064036	1	7	1065790

Wheels - Non-Projection Handrims and Attaching Hardware

Handrims, Composite

Model Number	Wheel Type	Handrim	Quantity	Attaching Points	Attaching Hardware
U2222	24" Composite Urethane	1063919	1	9	1065791

NOTE: Kits Contain Hardware for One Handrim