


Legrest Support Assembly

AS4A, AL4A/C, AHL4A/C

NOTE: Due to the possibility of various manufacturing dates for Tracer IV Wheelchairs built with a serial # starting with "03H" or "03I", replacement parts will need to be ordered thru an Invacare Parts Specialist at 1-800-333-6900.



Legrest Support Assembly

AS4A, AL4A/C, AHL4A/C

NOTE: Due to the possibility of various manufacturing dates for Tracer IV Wheelchairs built with a serial # starting with "03H" or "03I", replacement parts will need to be ordered thru an Invacare Parts Specialist at 1-800-333-6900.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		04526X019	Plate, Legrest Pad	1	1
2		1025104	Bracket, Legrest, Chrome	1	1
3		1022877-NLA	NLA - Tube, Support and Drive Link, Chrome - Right	1	1
3		1022878-NLA	NLA - Tube, Support and Drive Link, Chrome - Left	1	1
4		1022848	Tube, Slide, Articulating, Chrome - Right	1	1
4		1022849	Tube, Slide, Articulating, Chrome - Left	1	1
5		1027377	Washer (#20)	1	1
6		60264X001	Locknut (5/16-24)	1	1
7		11232X009	Bolt, Lug (5/16-24 x 1-3/8")	1	1
8		40363X000	Hex Nut (3/8"-16)	1	2
9		1022866-NLA	NLA - Screw, Socket Button Head (3/8-16 x 3/4")	1	1
10		1049527	Washer, Nylon (3/8 x 3/4 x 1/16")	1	2
11		1022846-NLA	NLA - Bar, Link, Black	1	1
12		1022874	Bumper	1	1
13		1030979	Screw, Phillips Pan Head Tap (#10-32 x 3/4")	1	1
14	A	1022871-NLA	NLA - Legrest Support Tube Assembly, Adult, Chrome - Right	1	1
14	A	1022872-NLA	NLA - Legrest Support Tube Assembly, Adult, Chrome - Left	1	1
14	B	1023092-NLA	NLA - Legrest Support Tube Assembly, Hemi, Chrome - Right	1	1
14	B	1023093-NLA	NLA - Legrest Support Tube Assembly, Hemi, Chrome - Left	1	1
15		1027097	Plug Button, Black (3/4")	1	1
16		1022854-NLA	NLA - Screw, Shoulder Slotted Head (1/2 x 1-1/8", 3/8-16)	1	1
17		1014807	Screw, Hex Head (5/16-24 x 3-1/4")	1	1
18		1091220	Spacer, Legrest	1	1
19		48823X008	Spring, Leaf	1	1
20		47328X003	Bushing, Pivot (5/16 x 7/16 x 7/8")	1	1
21		47330X003	Spacer (5/16 x 1/2 x 3/8")	1	1
22		10261X000	Nut, Acorn (5/16-24)	1	1
23		42326X004	Bar, Ratchet - Right/Left	1	1
24		04528X002	Bracket, Release, Black	1	1
25		1040060	Tip, Plastic, Black	1	1
26		1037125	Cam, Legrest	1	1
27		1037485	Spring, Release	1	1
28		1037126	Handle, Release	1	1
29		1013519	Washer, Nylon (5/16 x 11/16 x 1/16")	1	1
30		1040831	Ring, Snap	1	1
31		1107321	Nut, Low Cap (5/16"-24)	1	1
	C,D	1041799	Legrest Support Assembly, Articulating, Hemi, Chrome - Right	1	1

Legrest Support Assembly

AS4A, AL4A/C, AHL4A/C

NOTE: Due to the possibility of various manufacturing dates for Tracer IV Wheelchairs built with a serial # starting with "03H" or "03I", replacement parts will need to be ordered thru an Invacare Parts Specialist at 1-800-333-6900.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	C,D	1041800	Legrest Support Assembly, Articulating, Hemi, Chrome - Left	1	1
	C,E	1046329	Kit, Elevating Legrest Ratchet Hardware (Pair)	1	2
NOTE: A - Used on AS4A, AL4A, and AL4C. B - Used on AHL4A and AHL4C. C - Not shown. D - Includes items 1-31. E - Includes items 17-30 and instruction sheet.					