




Frames


	A4 Frame Orders
 NOTE: When replacing rigid frames for chairs made Before 8/95, a new seat frame is required.	Orders for an A4 frame must be handled as a custom order. When placing the order, the following information is required. 1. Tapered or "V" front 2. Width and Depth. 3. Front seat-to-floor. 4. Rear seat-to-floor. 5. Caster size. 6. Rear wheel size. 7. For fixed back specify degree of angle. 8. For fixed back specify height range. 9. Frame color.




Tapered Front & Adj. RSTF




"V" Front & Adj. RSTF




Tapered Front View



Fixed Back & Fixed RSTF
(Available w/ Tapered or "V" Front)



Fold Down Back & Fixed RSTF
(Available w/ Tapered or "V" Front)



"V" Front View

Frames

Frames - Chairs Built Before 11/26/96

Frame Front	Caster Size	Width	Front Seat to Floor	Adjustable RSTF	Fold Down Back Fixed RSTF	Fixed Back Fixed RSTF
Tapered	3"/5"	12-20"	17"	8881049328-NL A	8881060060-NL A	8881060069-NL A
Tapered	3"/5"	12-20"	19"	8881049329-NL A	8881060061-NL A	8881060070-NL A
Tapered	3"/5"	12-20"	21"	8881049330-NL A	8881060062-NL A	8881060071-NL A
"V"	3"/5"	12-20"	17"	8881056318-NL A	8881060063-NL A	8881060072-NL A
"V"	3"/5"	12-20"	19"	8881056319-NL A	8881060064-NL A	8881060073-NL A
"V"	3"/5"	12-20"	21"	8881056320-NL A	8881060065-NL A	8881060074-NL A
"V"	6"	12-20"	17"	8881057105-NL A	8881060066-NL A	8881060075-NL A
"V"	6"	12-20"	19"	8881057106-NL A	8881060067-NL A	8881060076-NL A
"V"	6"	12-20"	21"	8881057107-NL A	8881060068-NL A	8881060077-NL A

Frames - Chairs Built After 11/25/96

Adjustable RSTF	Fold Down Back Fixed RSTF	Fixed Back Fixed RSTF
8881070121	8881070122	8881070123