

Miscellaneous Accessories

NOTE: Due to the possibility of various manufacturing dates for Tracer IV Wheelchairs built with a serial # starting with "03H" or "03I", replacement parts will need to be ordered thru an Invacare Parts Specialist at 1-800-333-6900.

Single Pole Anti-Theft Device

Part Number	Description
1076646	1614 Single Pole Anti-Theft Device, Navy
1076647-NLA	1714 Single Pole Anti-Theft Device, Red

Anti-Tipper

Part Number	Description
1042412	Front (1326)
1036900	Rear, Swingaway and Non-Recliner Chairs (Not used on Chairs with One Arm Drives or Amputee Attachments.)
1036901-NLA	Rear, Chairs with One-Arm Drives and/or Amputee Attachments
1074064	Rear, Recliner Chairs (9758)
1061669-NLA	Rear, Fixed Front Frame

Carrying Straps

Part Number	Description
1045988-NLA	U05 Jet Black Lightweight Vinyl
8881091731-NLA See Note(s) A	All Other Colors

NOTE: A - Specify upholstery color, see page titled "Upholstery Colors".

Miscellaneous Accessories

NOTE: Due to the possibility of various manufacturing dates for Tracer IV Wheelchairs built with a serial # starting with "03H" or "03I", replacement parts will need to be ordered thru an Invacare Parts Specialist at 1-800-333-6900.

Crutch/Cane Carrier

Part Number	Description
1025467	Crutch and Cane Carrier

Non-Folding Device

Part Number	Description
1034205-NLA	Non-Folding Device (14" - 20" W)
1143048	Non-Folding Device (14" - 24" W)

Transfer Board

Part Number	Description
1034216	Transfer Board

OverHead Anti-Theft Device

Part Number	Description
1514	OverHead Anti-Theft Device Available on 16-18" Wide Chairs

Miscellaneous Accessories

NOTE: Due to the possibility of various manufacturing dates for Tracer IV Wheelchairs built with a serial # starting with "03H" or "03I", replacement parts will need to be ordered thru an Invacare Parts Specialist at 1-800-333-6900.

Seat/Chest Positioning Strap

Buckle Style	14-18" Wide Chair	20-22" Wide Chair	24" Wide Chair
Auto	1004310	1034625	1034622-NLA See Note(s) A
Airline	TAGTS1034751-NLA	1034752-NLA See Note(s) B	N/A
Hook and Loop	1034744	1034745	1034746
NOTE: A - For replacement contact The Aftermarket Group (TAG). B - 1034752 Airline Buckle Style is No Longer Available. Use Auto Style part number 1034625 as a replacement.			

Carry Pocket - Vinyl

Part Number	Description
8881047954-NLA See Note(s) A	14" Wide Chairs
8881047957 See Note(s) A	16 - 20" Wide Chairs
8881047960 See Note(s) A	22" Wide Chairs
NOTE: A - Specify upholstery color when ordering, see category titled, "Upholstery Colors."	

Miscellaneous Accessories

NOTE: Due to the possibility of various manufacturing dates for Tracer IV Wheelchairs built with a serial # starting with "03H" or "03I", replacement parts will need to be ordered thru an Invacare Parts Specialist at 1-800-333-6900.

Carrying Pocket - Sealed-On Back Upholstery

Upholstery Color	16-1/2" Wide Chair	18" Wide Chair	20" Wide Chair
U05 Jet Black Lightweight Vinyl	1045983-NLA	1045986-NLA	N/A
U60 Marina Blue Vinyl	1041608-NLA	1041610-NLA	1047154-NLA
U75 Dove Grey Vinyl	1041836-NLA	1041838-NLA	1047153-NLA
U85 Dark Blue Vinyl	1041614-NLA	1041616-NLA	1047150-NLA
U86 Ruby Vinyl	1041620-NLA	1041622-NLA	1047155-NLA
U88 Camel Vinyl	1041595-NLA	1041597-NLA	1047152-NLA
U92 Space/Sky Blue Vinyl	1041601-NLA	1041603-NLA	1047151-NLA

Solid Back Insert

Part Number	Description
8881041703-NLA See Note(s) A	16-1/2" Wide Chairs
8881034768-NLA See Note(s) A	18" Wide Chairs

NOTE: A - Specify upholstery color when ordering, see category titled, "Upholstery Colors."

Solid Seat Insert

Upholstery Type	16-1/2" Wide Chair	18" Wide Chair
Vinyl	8881041698-NLA See Note(s) A	8881034650 See Note(s) A
U05 Black Lightweight Vinyl	1046004-NLA	8881052548-NLA See Note(s) A

NOTE: A - Specify upholstery color when ordering, see category titled, "Upholstery Colors."