



Hideaway Undermount Wheel Lock (HUMW)



SCM_067A

Hideaway Undermount Wheel Lock (HUMW)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1093376-NLA	NLA - Wheel Lock Assembly, Undermount - Right	1	1
1		1093377-NLA	NLA - Wheel Lock Assembly, Undermount - Left	1	1
2		1006490	Tip, Plastic, Black	1	1
3		1068210-NLA	Kit, Wheel Lock Attaching Hardware (1" Tubing)	1	1
			Screw, Allen Head (1/4-20 x 1")	1	2
			Clamp, Half, Non-Threaded	1	1
			Clamp, Half, Threaded	1	1
4		1026096-NLA	NLA - Clamp, Wheel Lock, Threaded (1")	1	1
5		1026097-NLA	NLA - Clamp, Wheel Lock, Non-Threaded (1")	1	1
6		1027268	Screw, Allen Head (1/4-20 x 1")	1	2
NOTE: A - Includes Screws and Half Clamps.					