


Rigid Tubular Footplate Tapered Front Only



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1116379	Kit, Clamps Footrest	1	1
2			NLA - Clamp, Footrest	1	2
3			NLA - Screw, Socket Set Cone Point (1/4-20 x 1/2")	1	2
4		See Chart	#3 Rigid Tubular Footplate	1	1

#3 Rigid Tubular Footplate

Chair Width	Aluminum	Titanium
14"	1043272-NLA	1116359-NLA
15"	1043273-NLA	1116360-NLA
16"	1043274-NLA	1116361-NLA
17"	1043275-NLA	1116362-NLA
18"	1043276-NLA	1116363-NLA
19"	1043277-NLA	1116364-NLA
20"	1043278-NLA	1116365-NLA
21"	1043279-NLA	1116366-NLA
22"	1043280-NLA	1116367-NLA