



# Motor, Transformer and Pendant AC, DC - European / International / U.K.


**Motor AC and Pendant  
Before 5/1/98**



**Motor DC and Pendant  
From 5/1/98 thru 7/30/99**



**Motor,Transformer, and Pendant  
After 7/30/99**



| Item Number | Note | Part Number   | Description                                                                                 | Quantity Per |          |
|-------------|------|---------------|---------------------------------------------------------------------------------------------|--------------|----------|
|             |      |               |                                                                                             | Package      | Assembly |
| 1           | A    | 52000M336-NLA | NLA - Motor Assembly Kit, AC<br>(Includes: Motor, Bearings, Washers, and Instruction Sheet) | 1            | 1        |
| 1           | B    | 52000M340-NLA | NLA - Motor Assembly Kit, DC<br>(Includes: Motor, Bearings, Washers, and Instruction Sheet) | 1            | 1        |
| 1           | C    | 1091189-NLA   | NLA - Motor without Transformer                                                             | 1            | 1        |
| 2           | A    | 7001M098      | Pendant, AC Motor                                                                           | 1            | 1        |
| 2           | B    | 54000M755-NLA | NLA - Pendant, DC Motor                                                                     | 1            | 1        |
| 2           | B    | 54000M676-NLA | NLA - Pendant (May be used for DC Motor)                                                    | 1            | 1        |
| 2           | C    | 1092740-NLA   | NLA - Pendant, DC Motor                                                                     | 1            | 1        |
| 3           |      | 51000M562-NLA | NLA - Tie Wrap                                                                              | 1            | 1        |
| 4           | D    | 54000M802-NLA | NLA - Transformer, European or International Chairs                                         | 1            | 1        |
| 4           | D    | 54000M803-NLA | NLA - Transformer, U.K. Chairs                                                              | 1            | 1        |
| 5           |      | 1090519-NLA   | NLA - Motor, Transformer, and Pendant Kit - European/ International Chairs                  | 1            | 1        |
| 5           |      | 1090520-NLA   | NLA - Motor, Transformer, and Pendant Kit - U.K. Chairs                                     | 1            | 1        |

NOTE: A - Before 5/1/98  
 B - Between 5/2/98 through 7/30/99  
 C - After 7/31/99  
 D - Transformer is required only on chairs built after 7/30/99.