

Wheels without Handrims - 7/16" Axle (Chairs built After 6/23/99)

Wheel Style	Wheel Size	Tire Type	Single Wheel Assembly	Tire	Tube Insert	Bearing	Bearing Spacer
U1763BK	20x1-1/4"	Urethane	1038555 See Note(s) A, B	1037774	N/A	40117X024	1015756
1763BK	20 x 1-3/8"	Pneumatic	1034126 See Note(s) A	1017519	1017520	40117X024	1015756
4763BK	20 x 1-3/8"	Pneumatic w/ Flat Free Insert	1034127 See Note(s) A, B	1017519	1023322	40117X024	1015756
U1764BK	22 x 1-1/4"	Urethane	1032115 See Note(s) A	1028250	N/A	40117X024	1015756
1764BK	22 x 1-3/8"	Pneumatic	1027149 See Note(s) A	1018544	1018545	40117X024	1015756
4746BK	22 x 1-3/8"	Pneumatic w/ Flat Free Insert	1027150 See Note(s) B, A	1018544	1023323	40117X024	1015756
U2244BK 7 Spoke	24 x 1-1/4"	Urethane	1032117	1026904	N/A	40117X024	1015756
2244BK	24 x 1-3/8"	Pneumatic	1025123	40176X060	40176X061	40117X024	1015756
4244BK	24 x 1-3/8"	Pneumatic w/ Flat Free Insert	1025124	40176X060	1010906	40117X024	1015756
U2222 9-Spoke	24 x 1"	Urethane	1080378 See Note(s) B, F, E	1080400 See Note(s) F	N/A	40117X024	1015756
U2222C 9-Spoke	24 x 1"	Urethane Corded	1122726	1118209 See Note(s) G	N/A	40117X024	1015756

NOTE: A - When ordering 20" or 22" wheels, a wheel lock shoe spacer may also be required. See Wheel Locks.
B - Composite Wheels with Flat Free Inserts or Urethane Tires have a plug which caps an opening in the rim.
E - After 6/23/99
F - Non-Corded Tire replaced with a Corded Tire as of 11/16/09. Refer to Wheel Style U2222C part #1122726 for a wheel assembly or part #1118209 for a tire. For chairs manufactured before 11/16/09, it is recommended to order (2) each in quantity to match.
G - Corded Tires maybe more difficult to stretch over the rim. The complete wheel assembly that has the tire and wheel rim preassembled maybe advisable.

Description	Part Number
Rear Wheel Hardware	1047943