


Van Seat Frame and Mounting Hardware

9/1/02 - 7/13/03

Note: Seat post can be found on the base frame assembly page.



Van Seat Frame and Mounting Hardware

9/1/02 - 7/13/03

Note: Seat post can be found on the base frame assembly page.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	C	1113656-NLA	NLA - Frame, Van Seat, Black	1	1
1	D	1114349	Frame, Van Seat, Black	1	1
2		1026900	Washer (1/4 x 23/32 x 1/16")	1	6
3		1025713-NLA	NLA - Screw, Hex Head (1/4-20 x 1")	1	4
4		1114911-NLA	NLA - Knob, Adjustment, Black (5/16-18 x 1-1/2")	1	1
5	C	1114998-NLA	NLA - Knob Assembly, Black (3/8-16 x 9-1/4")	1	1
5	D	1114999-NLA	NLA - Knob Assembly, Black (3/8-16 x 10-1/4")	1	1
6		1088897-NLA	NLA - Screw, Hex Head w/ Patch (1/4-20 x 1-1/4")	1	2
7		1108410	Spacer (1/4 x 5/8 x 1/4")	1	2
8		1026099	Spacer, Coved (1/4 x 11/16 x 1/2")	1	4
9		1109472-NLA	NLA - Screw, Hex Head w/Patch (1/4-20 x 2-1/2")	1	2
10		1029501	Plug Button, Black (1")	1	2
11	A,C	1114794	Hardware, Van Seat Mounting w/ Plug Buttons	1	1
11	B,C	1114795-NLA	NLA - Hardware, Van Seat Mounting w/ Plug Buttons	1	1
11	A,D	1114792-NLA	NLA - Hardware, Wide Van Seat Mounting w/ Plug Buttons	1	1
11	B,D	1114793	Hardware, Wide Van Seat Mounting w/o Plug Buttons	1	1
12		1113691-NLA	NLA - Bracket, Center Post, Black	1	1

NOTE: A - Includes items 2-10
 B - Includes items 2-9
 C - 18" Wide Seat
 D - 20" Wide Seat