



Captain Seat Assembly (Before 11/1/01)


NOTE: Captain Seats incorporate an assisted recline of approximately 90 -114 degrees for extended backs and high back 22" wide models. High back 18" models are 90-170 degrees.



Extended Back 18" Wide




Extended Back, Forward Folding 18" Wide with Solid Seat Base



NOTE: High back 18" wide captain seat includes the actuator hardware. See seat frame and actuator hardware page.

High Back 18" Wide with Full Recline



High Back 22" Wide with Semi Recline

Captain Seat Assembly (Before 11/1/01)

NOTE: Captain Seats incorporate an assisted recline of approximately 90 -114 degrees for extended backs and high back 22" wide models. High back 18" models are 90-170 degrees.

Item Number	Backrest/Headrest Extensions and Captain Seat Assemblies	Grey Fabric Vinyl	All Grey Vinyl	Brunswick Buff All Tan Vinyl
1.	Backrest, Extended Back	1090374-NLA See Note(s) A	1090375-NLA See Note(s) A	1090376-NLA
1.	Backrest, Extended Back Forward Folding	1090374-NLA See Note(s) A	N/A	N/A
2.	Headrest, High Back (18 and 22" W)	1090379-NLA See Note(s) A	1090380-NLA See Note(s) A	N/A
3.	18" W Extended Back	1090371-NLA See Note(s) C	1090372-NLA See Note(s) C	1090373-NLA See Note(s) C
3.	18" W Extended Back Forward Folding	1090370-NLA See Note(s) C	N/A	N/A
3.	18" W High Back w/ Full Recline	1091068-NLA See Note(s) B, C	1091069-NLA See Note(s) C	N/A
3.	22" W High Back w/ Semi Recline	1098287-NLA	N/A	N/A

NOTE: A - 4-3/4" Backrest/Headrest Post Spacing
 B - High back 18" wide captain seat includes all actuator attaching hardware
 C - For captain seats marked NLA, use captain seats dated after 10/31/01

NOTE: When converting from an extended back captain seat to an 18" wide high back captain seat, the recliner wiring harness with sensor cable must be ordered as well.