


# PinDot ContourU Growth Links

GL01 Standard and GL03 Heavy Duty - Interfacing Hardware Package



# PinDot ContourU Growth Links

## GL01 Standard and GL03 Heavy Duty - Interfacing Hardware Package

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		GL01	Growth Links with Hardware, Standard (Pair)	2	1
1		GL03	Growth Links with Hardware, Heavy Duty (Pair)	2	1
2		3000634-NLA	NLA - Growth Link without Hardware, Standard - Left	1	2
2		3000635-NLA	NLA - Growth Link without Hardware, Standard - Right	1	2
2		3000627-NLA	NLA - Growth Link without Hardware, Heavy Duty - Left	1	2
2		3000628-NLA	NLA - Growth Link without Hardware, Heavy Duty - Right	1	2
3		1100023	Screw, Button Head w/ Patch (1/4-20 x 1/2")	1	10
4		1100039	Washer (1/4 x 9/16 x 1/16")	1	10
5		3000242-NLA	NLA - Backing Plate, 2 Position, 1" Spacing	1	4
	A,B	3000629-NLA	NLA - Growth Link Hardware	1	1
NOTE: A - Not Shown B - Includes items 3-5					