

## Rear Wheels

**NOTE: 5/8" Axle Wheels are used for fixed front frame models only.**

Wheel Style	Single Wheel Without Handrim	Tire	Tube/Insert	Bearing	Bearing Spacer	Flat Free Rim Plug
U2222 Composite Urethane 24 x 1" (9 Spoke) (5/8" Axle)	1080383-NLA  See Note(s) A	1080400  See Note(s) A	Not Needed	40117X006	1019695-NLA	Not Needed
U2222C Composite Urethane 24 x 1" (9 Spoke) Low Profile (5/8" Axle)	1121043	1118209  See Note(s) B	Not Needed	40117X006	1019695-NLA	Not Needed
U2222 Composite Urethane 24 x 1" (9 Spoke) Low Profile (7/16" Axle)	1080378  See Note(s) C	1080400  See Note(s) C	Not Needed	40117X024	1015756	Not Needed
U2222C Composite Urethane 24 x 1" (9 Spoke) Low Profile (7/16" Axle)	1122726	1118209  See Note(s) B	Not Needed	40117X024	1015756	Not Needed
2244BK Composite Pneumatic 24 x 1-3/8" (5/8" Axle)	1034128	40176X060	40176X061	40117X006	1019695-NLA	Not Needed
4244BK Composite Pneumatic w/Flat Free 24 x 1-3/8" (5/8" Axle)	1034129	40176X060	Not Needed	40117X006	1019695-NLA	Not Needed
U1763BK Composite Urethane 20 x 1-1/4" (7/16" Axle)	1038555	1037774	Not Needed	40117X024	1015756	Not Needed
U1764BK Composite Urethane 22 x 1-3/8" (7/16" Axle)	1032115	1028250	Not Needed	40117X024	1015756	Not Needed

## Rear Wheels

**NOTE: 5/8" Axle Wheels are used for fixed front frame models only.**

Wheel Style	Single Wheel Without Handrim	Tire	Tube/Insert	Bearing	Bearing Spacer	Flat Free Rim Plug
2244BK Composite Pneumatic 24 x 1-3/8" (7/16" Axle)	1025123	40176X060	40176X061	40117X024	1015756	Not Needed
4244BK Composite Pneumatic w/Flat Free 24 x 1-3/8" (7/16" Axle)	1025124	40176X060	1010906	40117X024	1015756	1030936
U2244BK Composite Urethane (7 Spoke) 24 x 1-1/4" (7/16" Axle)	1032117	1026904	Not Needed	40117X024	1015756	Not Needed
<p>NOTE: A - Non-Corded Tire replaced with a Corded Tire as of 11/16/09. Refer to Wheel Style U2222C part #1121043 for a wheel assembly or part #1118209 for a tire. For chairs manufactured before 11/16/09, it is recommended to order (2) each in quantity to match</p> <p>B - Corded Tires maybe more difficult to stretch over the wheel rim. The complete wheel assembly that has the tire and wheel rim preassembled maybe advisable.</p> <p>C - Non-Corded Tire replaced with a Corded Tire as of 11/16/09. Refer to Wheel Style U2222C part #1122726 for a wheel assembly or part #1118209 for a tire. For chairs manufactured before 11/16/09, it is recommended to order (2) each in quantity to match</p>						