

## Sling Type and Fabric

NOTE : All Slings are sold as individual units unless otherwise stated.

Model Number	Description	Fabric Type
9040NLA	Wrap Around with out Headrest and with Commode Opening	Navy Blue Mesh
❖⊕ 9042	One Piece without Headrest and Commode Opening	Blue Polyester
9042PNLA	One Piece without Headrest and Commode Opening	Blue Polyester
⊕ 9042XWNLA	One Piece without Headrest and Commode Opening	White Polyester
❖⊕ 9043	One Piece without Headrest and with Commode Opening	Blue Polyester
9043PNLA	One Piece without Headrest and with Commode Opening	Blue Polyester
⊕ 9043XWNLA	One Piece without Headrest and with Commode Opening	Blue Polyester
⊕ 9044NLA	One Piece with Headrest and with out Commode Opening	Blue Polyester
⊕ 9045NLA	One Piece with Headrest and Commode Opening	Blue Polyester
❖⊕ 9046	One Piece without Headrest and Commode Opening	Navy Blue Mesh
⊕ 9046PNLA	One Piece without Headrest and Commode Opening	Navy Blue Mesh
9046XWNLA	One Piece without Headrest and Commode Opening	White Mesh
❖⊕ 9047	One Piece without Headrest and with Commode Opening	Nylon Mesh
⊕ 9047PNLA	One Piece without Headrest and with Commode Opening	Navy Blue Mesh
⊕ 9047XWNLA	One Piece without Headrest and with Commode Opening	Blue Polyester
⊕ 9048NLA	One Piece with Headrest and with out Commode Opening	Blue Polyester
⊕ 9049NLA	One Piece with Headrest and Commode Opening	Blue Polyester
9742NLA	Wrap Around with out Headrest and Commode Opening	Teal Mesh
9746NLA	One Piece with Headrest and with out Commode Opening	Teal Mesh
9747NLA	One Piece with Headrest and Commode Opening	Teal Mesh
R100P	Divided Leg Sling - Petite	Navy Blue Polyester
R100	Divided Leg Sling - Medium	Navy Blue Polyester
R101	Divided Leg Sling - Large	Navy Blue Polyester
R102	Divided Leg Sling - X-Large	Navy Blue Polyester
R110	Full Body Sling without Commode Opening - Medium	Navy Blue Mesh Poly.
R111	Full Body Sling without Commode Opening - Large	Navy Blue Mesh Poly.
R112	Full Body Sling without Commode Opening - Medium	Navy Blue Polyester
R113	Full Body Sling without Commode Opening - Large	Navy Blue Polyester
R114	Full Body Sling with Commode Opening - Medium	Navy Blue Mesh Poly.
R115	Full Body Sling with Commode Opening - Large	Navy Blue Mesh Poly.
R116	Full Body Sling with Commode Opening - X-Large	Navy Blue Mesh Poly.
R117	Full Body Sling without Commode Opening - X-Large	Navy Blue Mesh Poly.
R120NLA	Toileting Sling - Medium	Navy Blue Polyester
R121	Toileting Sling - Large	Navy Blue Polyester
R122NLA	Toileting Sling - X-Large	Navy Blue Polyester
R130	Stand-up Sling (RPS350 Stand-up Lift Only)	Navy Blue Polyester
R134	Transport Sling (RPS350 Stand-up Lift Only)	Navy Blue Polyester
R140	Full Body Heavy Duty Sling with out Commode Opening	Navy Blue Mesh Poly.
R141	Full Body Heavy Duty Sling with Commode Opening	Navy Blue Mesh Poly.

❖ These products are used in connection with CareGuard™ by Invacare®.

⊕ Sling manufactured before 8/25/97 have been made from a different sling material.

## Sling Measurements

NOTE : All Slings are sold as individual units unless otherwise stated.

⊛ Model Number	Overall Length	Overall Width	Headrest	Commode Opening To Top of Back	Commode Opening To Front of Seat	Commode Opening	
						Width	Depth
⊛ 9040NLA	36"	31-1/2"	—	23-1/4"	—	7"	8"
❖⊛ 9042	40-1/2"	29"	—	—	—	—	—
9042PNLA	32-1/2"	26-1/2"	—	—	—	—	—
⊛9042XWNLA	46-1/2"	41"	—	—	—	—	—
❖⊛ 9043	40-1/2"	29"	—	25-1/2"	4-1/2"	6-1/2"	11"
9043PNLA	32-1/2"	26-1/2"	—	17"	4-1/2"	6-1/2"	11"
⊛ 9043XWNLA	46-1/2"	41-1/2"	—	29"	6-1/2"	6-1/2"	11"
⊛ 9044NLA	54-1/2"	29"	14"	—	—	—	—
⊛ 9045NLA	54-1/2"	29"	14"	25-1/2"	4-1/2"	6-1/2"	11"
❖⊛ 9046	40-1/2"	29"	—	—	—	—	—
⊛ 9046PNLA	32-1/2"	26-1/2"	—	—	—	—	—
9046XWNLA	46-1/2"	41"	—	—	—	—	—
❖★⊛ 9047	40-1/2"	29"	—	25-1/2"	4-1/2"	6-1/2"	11"
★⊛ 9047PNLA	32-1/2"	26-1/2"	—	17"	4-1/2"	5-1/2"	10"
★⊛ 9047XWNLA	46-1/2"	41-1/2"	—	29"	6-1/2"	6-1/2"	11"
⊛ 9048NLA	54-1/2"	29"	14"	—	—	—	—
⊛ 9049NLA	54-1/2"	29"	14"	25-1/2"	4-1/2"	6-1/2"	11"
9742NLA	52"	41"	—	—	—	—	—
9746NLA	50"	40"	—	—	—	—	—
9747NLA	50"	40"	—	—	—	6"	12"
R100P	50"	33"	—	—	—	—	—
R100	53-1/2"	37-1/2"	—	—	—	—	—
R101	58-1/2"	41"	—	—	—	—	—
R102	63-1/2"	44-1/2"	—	—	—	—	—
R110	53-1/4"	43-3/4"	—	—	—	—	—
R111	58-1/4"	47-3/4"	—	—	—	—	—
R112	53-1/4"	43-3/4"	—	—	—	—	—
R113	58-1/4"	47-3/4"	—	—	—	—	—
R114	53-1/4"	43-3/4"	—	34-3/4"	8"	6"	11"
R115	58-1/4"	47-3/4"	—	38"	7"	8"	13-1/2"
R116	63-1/4"	47-3/4"	—	43"	7"	7"	13"
R117	63-1/2"	47-3/4"	—	—	—	—	—
R120NLA	33"	41-1/4"	—	—	—	—	—
R121	37"	44-3/4"	—	—	—	—	—
R122NLA	42"	48-1/4"	—	—	—	—	—
R130	15"	39"	—	—	—	—	—
R134	39"	29"	—	—	—	—	—
R140	55"	45"	—	—	—	—	—
R141	55"	45"	—	41-1/2"	13-1/2"	7"	13"

⊛ Sling manufactured before 8/25/97 have been made from a different sling material.

❖ These products are used in connection with CareGuard™ by Invacare®.

★ Commode Opening may vary in width and/or depth 1/2" to 1" in size.