


Seat Frame

THE PRODUCTION OF THE HTR IS CURRENTLY SUSPENDED AND SOME SERVICE PARTS MAY NOT BE AVAILABLE TO ORDER.



Seat Frame

THE PRODUCTION OF THE HTR IS CURRENTLY SUSPENDED AND SOME SERVICE PARTS MAY NOT BE AVAILABLE TO ORDER.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		460068-NLA	NLA - Pan, Seat, with Hook and Loop Fasteners (16" and 18" W)	1	1
2		1108191-NLA	NLA - Screw, Flat Head (1/4-10 x 1")	1	4
2		0288706-NLA	NLA - Screw, Socket, Flat Head (1/4-20 x 1")	1	4
3		1108106	Washer, Nylon (5/16 x 3/4 x 1/8")	1	8
4		1025736	Locknut (1/4-20)	1	4
5		450088-25-NLA	Seat, Contoured, Mauve (16")	1	1
5		450089-25-NLA	Seat, Contoured, Mauve (18")	1	1
5		450088-26-NLA	Seat, Contoured, Grotto (16")	1	1
5		450089-26-NLA	Seat, Contoured, Grotto (18")	1	1
5		450088-27-NLA	Seat, Contoured, Blueberry (16")	1	1
5		450089-27-NLA	NLA - Seat, Contoured, Blueberry (18")	1	1
6		1108194-NLA	NLA - Locknut (#10-32)	1	2
7		1108134	Knob, Black	1	2
8		TR46B	Sleeve, Plastic	1	1
9		1013560-NLA	NLA - Screw, Hex Head (5/16-18 x 1-1/2")	1	4
10		440544	Slide, Seat - Left	1	1
11		1040292-NLA	NLA - Washer, Steel Flat (1/4 x 5/8 x 1/16")	1	4
12		1108192-NLA	Screw, Hex Head (1/4-20 x 1")	1	4
13		440545	Slide, Seat - Right	1	1
14		420131-1-NLA	Adjuster, Leg (16")	1	1
14		420132-1	Adjuster, Leg (18")	1	1
15		16308X000	Washer, Flat (5/16 x 11/16 x 1/16")	1	4
16		1108122-NLA	NLA - Locknut (5/16-18)	1	4
17		60290X002	Locknut (5/16-18)	1	2
18		1108114-NLA	NLA - Washer, Nylon (5/16 x 1/2 x 1/16")	1	2
19		1013560-NLA	NLA - Screw, Hex Head (5/16-18 x 1-1/2")	1	1
20		1108188-NLA	Spring, Lift Pin	1	2
21		1132722	Pin, Release	1	2
22	J	350025-NLA	NLA - Strut, Gas	1	1
22		1047705	Cylinder, Lock, Tilt - Left	1	1
23		350039-NLA	NLA - Cable, Tilt with Handle	1	1
24		TR432	Activator, Strut	1	1
25		1108155-NLA	Clip, Tinnerman, U-Type	1	1
26		16325X000-NLA	NLA - Washer (3/8 x 13/16 x 1/16")	1	1
27		B444-NLA	NLA - Spacer (5/16 x 7/16 x 3/4")	1	1
28		1108121-NLA	NLA - Screw, Hex Head (5/16-18 x 2")	1	1
29	G	440492-1-NLA	Frame, Seat, Chrome (16")	1	1
29	H	1129532-NLA	Frame, Seat, Grey (16")	1	1
29	G	440493-1-NLA	Frame, Seat, Chrome (18")	1	1
29	H	1129533-NLA	Frame, Seat, Grey (18")	1	1
30		1071237	Plug Button, Square, Black (1")	1	2
31		1088128	Guard, Finger, Tilt	1	2
	A,B	440569	Strut Assembly	1	1

Seat Frame

THE PRODUCTION OF THE HTR IS CURRENTLY SUSPENDED AND SOME SERVICE PARTS MAY NOT BE AVAILABLE TO ORDER.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A,C	440575	Seat Base Assembly	1	1
	A,D	440580	Leg Adjuster Assembly (16")	1	1
	A,D	440581	Leg Adjuster Assembly (18")	1	1
	A,E	440564	Kit, Seat Frame Hardware	1	1
			NLA - Screw, Hex Head (5/16-18 x 1-1/2")	1	4
			Plug Button, Square, Black (1")	1	2
			Washer, Nylon (5/16 x 3/4 x 1/8")	1	8
			NLA - Locknut (5/16-18)	1	4
			Knob, Black	1	2
			Spring, Lift Pin	1	2
			Washer, Flat (5/16 x 11/16 x 1/16")	1	4
			Slide, Seat - Left	1	1
			Slide, Seat - Right	1	1
			Pin, Release	1	2
			Sleeve, Plastic	1	1
	A,F	440574	Kit, Contoured Seat Hardware	1	1
			NLA - Washer, Steel Flat (1/4 x 5/8 x 1/16")	1	4
			Screw, Hex Head (1/4-20 x 1")	1	4
<p>NOTE: A - Not Shown. B - Includes items 17-19, 22-25 and 27-28. C - Includes items 1-4. D - Includes items 6-7, 14 and 20-21. E - Kit Includes Items 7-10, 13, 15-16, 20-21 and 30. F - Kit Includes Items 11-12. G - Before 1/24/05. For use on chairs with "Chrome Frame Finish." H - After 1/23/05. For use on chairs with "Grey Frame Finish." J - No longer available as a service part, use 1047705 as the replacement part for this item.</p>					