


Omit Recline Tube

THE PRODUCTION OF THE HTR IS CURRENTLY SUSPENDED AND SOME SERVICE PARTS MAY NOT BE AVAILABLE TO ORDER.



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	440583-NLA	NLA - Omit Recline Tube Assembly - Chrome	1	1
1	B	1129599-NLA	Omit Recline Tube Assembly - Grey	1	1
2	A	440578-NLA	NLA - Tube, Omit Recline, Chrome	1	1
2	B	1129541-NLA	Tube, Omit Recline, Grey	1	1
3		1108121-NLA	NLA - Screw, Hex Head (5/16-18 x 2")	1	2
4		60290X002	Locknut (5/16-18)	1	2

NOTE: A - Before 1/24/05. For use on chairs with "Chrome Frame Finish."
 B - After 1/23/05. For use on chairs with "Grey Frame Finish."