


Forks & Hardware

THE PRODUCTION OF THE HTR IS CURRENTLY SUSPENDED AND SOME SERVICE PARTS MAY NOT BE AVAILABLE TO ORDER.



*6" Caster and Frame are
Shown for Reference only*

Forks & Hardware

THE PRODUCTION OF THE HTR IS CURRENTLY SUSPENDED AND SOME SERVICE PARTS MAY NOT BE AVAILABLE TO ORDER.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1096414-NLA	NLA - Cap, End, Black Plastic	1	1
2		1045608	Locknut (7/16-20)	1	1
3		1006049	Bearing (7/16")	1	2
4		1046811	Fork, 8" Caster (7/16" Stem)	1	1
4		1019972	Fork, 6" Caster (7/16" Stem)	1	1
5	B	430053-NLA	NLA - Fork Assembly (8" Caster)	1	1
5	B	1088016	Fork Assembly (6" Caster)	1	1

NOTE: A - As of 8/22/00, End Cap changed from Nickel Plated to Black Plastic.
 B - Includes items 1-4.