




Captain Seat Assembly (Before 9/1/02)



Extended Back 18" Wide



Extended Back, Forward Folding 18" Wide with Solid Seat Base



NOTE: High back 18" wide captain seat includes the actuator hardware.

High Back 18" Wide with Full Recline

Captain Seat Assembly (Before 9/1/02)

Item #1 - Extended Back Backrest

Extended Back Backrest (E.B.)	Grey Fabric/Vinyl	All Grey Vinyl	Brunswick Buff All Tan Vinyl
E.B. Before 8/4/99 and Forward Folding Before 5/12/00	1078997-NLA See Note(s) A	N/A	N/A
E.B. After 8/3/99 and Forward Folding After 5/11/00	1090374-NLA See Note(s) B	N/A	N/A
E.B. Before 7/22/99	N/A	1081285-NLA See Note(s) A	1086109-NLA See Note(s) A
E.B. After 7/21/99	N/A	1090375-NLA See Note(s) B	1090376-NLA See Note(s) B
NOTE: A - 9" Backrest Post Spacing. B - 4-3/4" Backrest/Headrest Post Spacing.			

Item #2 - High Back Headrest

Headrest	Grey Fabric/Vinyl	All Grey Vinyl	Brunswick Buff All Tan Vinyl
Headrest, High Back Before 8/16/99	1071576-NLA See Note(s) C	N/A	N/A
Headrest, High Back After 8/15/99	1090379-NLA See Note(s) B	N/A	N/A
Headrest, High Back	N/A	1090380-NLA See Note(s) C	N/A
NOTE: B - 4-3/4" Backrest/Headrest Post Spacing. C - 4-23/32" Headrest Post Spacing.			

Captain Seat Assembly (Before 9/1/02)

Item #3 - 18" Wide Extended Back Captain Seat

Captain Seat - Extended Back	Grey Fabric/Vinyl	All Grey Vinyl	Brunswick Buff All Tan Vinyl
Before 8/4/99	1078965-NLA See Note(s) D	N/A	N/A
After 8/3/99	1090371-NLA See Note(s) E, D	N/A	N/A
Before 7/22/99	N/A	1081284-NLA See Note(s) D	1086108-NLA See Note(s) D
After 7/21/99	N/A	1090372-NLA See Note(s) E, D	1090373-NLA See Note(s) E, D
and Forward Folding Before 5/12/00	1085104-NLA See Note(s) D	N/A	N/A
and Forward Folding After 5/11/00	1090370-NLA See Note(s) D, E	N/A	N/A

NOTE: D - For Captain Seats marked NLA, use Captain Seats dated after 8/31/02 as replacements.
E - Backrest/Headrest post spacing changed. When replacing seat assemblies made before the date specified, a new style backrest or headrest must be ordered as well to match.

NOTE: When converting between high back and extended back captain seats, the wiring harness will also need to be converted. See page titled "Battery Charger and Harnesses".

Item #3 - 18" Wide High Back Captain Seat

Captain Seat - High Back	Grey Fabric/Vinyl	All Grey Vinyl	Brunswick Buff All Tan Vinyl
18" Wide High Back Before 8/16/99	1070456-NLA See Note(s) F, D	N/A	N/A
18" Wide High Back After 8/15/99	1091068-NLA See Note(s) E, F, D	N/A	N/A
18" Wide High Back	N/A	1091069-NLA See Note(s) F, D	N/A

NOTE: D - For Captain Seats marked NLA, use Captain Seats dated after 8/31/02 as replacements.
E - Backrest/Headrest post spacing changed. When replacing seat assemblies made before the date specified, a new style backrest or headrest must be ordered as well to match.
F - High back 18" wide captain seats includes all actuator attaching hardware (LESS) the limit switch sensor harness #1074983. See page titled "Captain Seat Frame and Actuator Hardware".

NOTE: When converting between high back and extended back captain seats, the wiring harness will also need to be converted. See page titled "Battery Charger and Harnesses".