


Axle Tube Position Assembly


For Suspension Chairs, it is recommended that Axle Tube Position Assemblies be replaced in pairs. First, select a Rear Suspension Assembly. Second, select a Positioning Tube.

When ordering, please state if the Rear Suspension Assembly and Positioning Tubes are to be factory assembled.


1. Rear Suspension Assembly



2. Positioning Tube for Suspension Chairs



3. Positioning Tube for Non-Suspension Chairs



Axle Tube Position Assembly

For Suspension Chairs, it is recommended that Axle Tube Position Assemblies be replaced in pairs. First, select a Rear Suspension Assembly. Second, select a Positioning Tube.

When ordering, please state if the Rear Suspension Assembly and Positioning Tubes are to be factory assembled.

#1 Rear Suspension Assembly

User's Weight Range	Part Number	Spring Color
0 lbs to 100 lbs	1093150-NLA	Blue/Blue
101 lbs to 140 lbs	1093151-NLA	Blue/Black
141 lbs to 180 lbs	1093152-NLA	Black/Black
181 lbs to 220 lbs	1093153-NLA	Black/Purple
221 lbs to 250 lbs	1093154-NLA	Purple/Purple

Axle Tube Position Assembly

For Suspension Chairs, it is recommended that Axle Tube Position Assemblies be replaced in pairs. First, select a Rear Suspension Assembly. Second, select a Positioning Tube.

When ordering, please state if the Rear Suspension Assembly and Positioning Tubes are to be factory assembled.

#2 Positioning Axle Tubes - Chairs without Amputee Attachments - Suspension Chairs

Wheel Size/ Chair Depth	17" Front STF	18" Front STF	19" Front STF	20" Front STF	21" Front STF
22" / 14" Deep	1098018-NLA	1098018-NLA	1098018-NLA	1098016-NLA	N/A
22" / 15" Deep	1098020-NLA	1098020-NLA	1098020-NLA	1098018-NLA	N/A
22" / 16" Deep	1098022-NLA	1098022-NLA	1098022-NLA	1098020-NLA	N/A
22"/ 17-20" Deep	1098024-NLA	1098024-NLA	1098024-NLA	1098022-NLA	N/A
24" / 14" Deep	1098020-NLA	1098018-NLA	1098018-NLA	1098018-NLA	1098016-NLA
24" / 15" Deep	1098022-NLA	1098020-NLA	1098020-NLA	1098020-NLA	1098018-NLA
24" / 16" Deep	1098024-NLA	1098022-NLA	1098022-NLA	1098022-NLA	1098020-NLA
24"/ 17-20" Deep	1098026-NLA	1098024-NLA	1098024-NLA	1098024-NLA	1098022-NLA
26" / 14" Deep	N/A	1098020-NLA	1098018-NLA	1098018-NLA	1098018-NLA
26" / 15" Deep	N/A	1098022-NLA	1098020-NLA	1098020-NLA	1098020-NLA
26" / 16" Deep	N/A	1098024-NLA	1098022-NLA	1098022-NLA	1098022-NLA
26"/ 17-20" Deep	N/A	1098026-NLA	1098024-NLA	1098024-NLA	1098024-NLA

Axle Tube Position Assembly

For Suspension Chairs, it is recommended that Axle Tube Position Assemblies be replaced in pairs. First, select a Rear Suspension Assembly. Second, select a Positioning Tube.

When ordering, please state if the Rear Suspension Assembly and Positioning Tubes are to be factory assembled.

#2 Positioning Axle Tubes - Chairs with Amputee Attachments - Suspension Chairs

Wheel Size/ Chair Depth	17" Front STF	18" Front STF	19" Front STF	20" Front STF	21" Front STF
22" / 14" Deep	1098026-NLA	1098026-NLA	1098026-NLA	1098024-NLA	N/A
22" / 15" Deep	1098028-NLA	1098028-NLA	1098028-NLA	1098026-NLA	N/A
22" / 16" Deep	1098030-NLA	1098030-NLA	1098030-NLA	1098028-NLA	N/A
22"/ 17-20" Deep	1098032-NLA	1098032-NLA	1098032-NLA	1098030-NLA	N/A
24" / 14" Deep	1098028-NLA	1098026-NLA	1098026-NLA	1098026-NLA	1098024-NLA
24" / 15" Deep	1098030-NLA	1098028-NLA	1098028-NLA	1098028-NLA	1098026-NLA
24" / 16" Deep	1098032-NLA	1098030-NLA	1098030-NLA	1098030-NLA	1098028-NLA
24"/ 17-20" Deep	1098034-NLA	1098032-NLA	1098032-NLA	1098032-NLA	1098030-NLA
26" / 14" Deep	N/A	1098028-NLA	1098026-NLA	1098026-NLA	1098026-NLA
26" / 15" Deep	N/A	1098030-NLA	1098030-NLA	1098028-NLA	1098028-NLA
26" / 16" Deep	N/A	1098032-NLA	1098032-NLA	1098030-NLA	1098030-NLA
26"/ 17-20" Deep	N/A	1098034-NLA	1098034-NLA	1098032-NLA	1098032-NLA

Axle Tube Position Assembly

For Suspension Chairs, it is recommended that Axle Tube Position Assemblies be replaced in pairs. First, select a Rear Suspension Assembly. Second, select a Positioning Tube.

When ordering, please state if the Rear Suspension Assembly and Positioning Tubes are to be factory assembled.

#3 Positioning Axle Tubes - Chairs without Amputee Attachments - Non-Suspension Chairs

Wheel Size/ Chair Depth	17" Front STF	18" Front STF	19" Front STF	20" Front STF	21" Front STF
22" / 14" Deep	1098088-NLA	1098088-NLA	1098086-NLA	1098086-NLA	N/A
22" / 15" Deep	1098088-NLA	1098088-NLA	1098086-NLA	1098086-NLA	N/A
22" / 16" Deep	1098088-NLA	1098088-NLA	1098086-NLA	1098086-NLA	N/A
22"/ 17-20" Deep	1098088-NLA	1098088-NLA	1098086-NLA	1098086-NLA	N/A
24" / 14" Deep	1098088-NLA	1098088-NLA	1098088-NLA	1098086-NLA	1098086-NLA
24" / 15" Deep	1098088-NLA	1098088-NLA	1098088-NLA	1098086-NLA	1098086-NLA
24" / 16" Deep	1098088-NLA	1098088-NLA	1098088-NLA	1098086-NLA	1098086-NLA
24"/ 17-20" Deep	1098088-NLA	1098088-NLA	1098088-NLA	1098086-NLA	1098086-NLA
26" / 14" Deep	N/A	1098088-NLA	1098088-NLA	1098088-NLA	1098086-NLA
26" / 15" Deep	N/A	1098088-NLA	1098088-NLA	1098088-NLA	1098086-NLA
26" / 16" Deep	N/A	1098088-NLA	1098088-NLA	1098088-NLA	1098086-NLA
26"/ 17-20" Deep	N/A	1098088-NLA	1098088-NLA	1098088-NLA	1098086-NLA

Axle Tube Position Assembly

For Suspension Chairs, it is recommended that Axle Tube Position Assemblies be replaced in pairs. First, select a Rear Suspension Assembly. Second, select a Positioning Tube.

When ordering, please state if the Rear Suspension Assembly and Positioning Tubes are to be factory assembled.

#3 Positioning Axle Tubes - Chairs with Amputee Attachments - Non-Suspension Chairs

Wheel Size/ Chair Depth	17" Front STF	18" Front STF	19" Front STF	20" Front STF	21" Front STF
22" / 14" Deep	1098096-NLA	1098096-NLA	1098094-NLA	1098094-NLA	N/A
22" / 15" Deep	1098096-NLA	1098096-NLA	1098094-NLA	1098094-NLA	N/A
22" / 16" Deep	1098096-NLA	1098096-NLA	1098094-NLA	1098094-NLA	N/A
22"/ 17-20" Deep	1098096-NLA	1098096-NLA	1098094-NLA	1098094-NLA	N/A
24" / 14" Deep	1098096-NLA	1098096-NLA	1098096-NLA	1098094-NLA	1098094-NLA
24" / 15" Deep	1098096-NLA	1098096-NLA	1098096-NLA	1098094-NLA	1098094-NLA
24" / 16" Deep	1098096-NLA	1098096-NLA	1098096-NLA	1098094-NLA	1098094-NLA
24"/ 17-20" Deep	1098026-NLA	1098096-NLA	1098096-NLA	1098094-NLA	1098094-NLA
26" / 14" Deep	N/A	1098096-NLA	1098096-NLA	1098094-NLA	1098094-NLA
26" / 15" Deep	N/A	1098096-NLA	1098096-NLA	1098094-NLA	1098094-NLA
26" / 16" Deep	N/A	1098096-NLA	1098096-NLA	1098094-NLA	1098094-NLA
26"/ 17-20" Deep	N/A	1098096-NLA	1098096-NLA	1098094-NLA	1098024-NLA