


# One Piece Elevating Footboard

THE PRODUCTION OF THE HTR IS CURRENTLY SUSPENDED AND SOME SERVICE PARTS MAY NOT BE AVAILABLE TO ORDER.



# One Piece Elevating Footboard

**THE PRODUCTION OF THE HTR IS CURRENTLY SUSPENDED AND SOME SERVICE PARTS MAY NOT BE AVAILABLE TO ORDER.**

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	D	420145-1-NLA	Tube, Frame, Footboard, Chrome (16")	1	1
1	E	1129510-NLA	Footboard, Frame, Tube, Grey (16")	1	1
1	D	420146-1	Tube, Frame, Footboard, Chrome (18")	1	1
1	E	1129511	Footboard, Frame, Tube, Grey (18")	1	1
2		420148	Footboard - ABS (16")	1	1
2		420149-NLA	NLA - Footboard - ABS (18")	1	1
3	D	420122-1-NLA	Footrest Lower, Support Tube, Chrome (16" W)	1	1
3	E	1129525-NLA	Footrest Lower, Support Tube, Grey (16" W)	1	1
3	D	420123-1-NLA	Footrest Lower, Support Tube, Chrome (18" W)	1	1
3	E	1129526-NLA	Footrest Lower, Support Tube, Grey (18" W)	1	1
4		1108128	Screw, Socket Button Head (1/4-20 x 1-1/2")	1	4
5		1025195	Locknut (1/4-20)	1	2
6		420131-1-NLA	Adjuster, Leg (16")	1	1
6		420132-1	Adjuster, Leg (18")	1	1
7		1108134	Knob, Black	1	2
8		1132724-NLA	NLA - Pin, Lift	1	1
9		1132722	Pin, Release	1	1
10		1108188-NLA	Spring, Lift Pin	1	2
11		1108194-NLA	NLA - Locknut (#10-32)	1	2
12		1108117	Spacer, Ribbed (25/32 x 1/4 x 7/8")	1	2
13	A,D	420182	Hanger Assembly, Telescoping, Chrome (16")	1	1
13	A,E	1131375	Hanger Assembly, Telescoping, Grey (16")	1	1
13	A,D	420183	Hanger Assembly, Telescoping, Chrome (18")	1	1
13	A,E	1129595	Hanger Assembly, Telescoping, Grey (18")	1	1
13	D	420112-1-NLA	Hanger, Telescoping, Chrome (16")	1	1
13	E	1129523-NLA	Hanger, Telescoping, Grey (16")	1	1
13	D	420113-1-NLA	Hanger, Telescoping, Chrome (18")	1	1
13	E	1129524-NLA	Hanger, Telescoping, Grey (18")	1	1
14	C	1108160-NLA	NLA - Block, Glide	1	2
14		1108498-NLA	NLA - Plug Button, Square, Black (1-1/4")	1	2
15		420147-NLA	NLA - Handle, Chrome	1	1
16	D	420142-1-NLA	Tube, Lower Support, Chrome	1	2
16	E	1129527-NLA	NLA - Tube, Lower Support, Grey	1	2
17		1008882	Locknut (1/4-20)	1	2
18		350046	Cylinder, Lock	1	1
19		1108192-NLA	Screw, Hex Head (1/4-20 x 1")	1	1
20		1108123-NLA	NLA - Screw, Hex Head (5/16-18 x 1-3/4")	1	1
21		60290X002	Locknut (5/16-18)	1	1
22		60275X000	Locknut (1/4-20)	1	1
23		350038-NLA	NLA - Cable	1	1
24		1035056-NLA	NLA - Screw, Flat Head (#10 x 1")	1	4
25		1010398	Washer, Nylon (1/4 x 1/2 x 1/16")	1	8
26		1108108-NLA	NLA - Button, Lock	1	2
27		420174-25-NLA	Pad, Calf, Lower, Mauve (16" W)	1	1

# One Piece Elevating Footboard

**THE PRODUCTION OF THE HTR IS CURRENTLY SUSPENDED AND SOME SERVICE PARTS MAY NOT BE AVAILABLE TO ORDER.**

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
27		420174-26-NLA	Pad, Calf, Lower, Grotto (16" W)	1	1
27		420174-27-NLA	Pad, Calf, Lower, Blueberry (16" W)	1	1
27		420175-25-NLA	Pad, Calf, Lower, Mauve (18" W)	1	1
27		420175-26-NLA	Pad, Calf, Lower, Grotto (18" W)	1	1
27		420175-27	Pad, Calf, Lower, Blueberry (18" W)	1	1
28	F	1121008-NLA	NLA - Screw (#10-24 x 7/8")	1	4
29	B,D	1072296	Footboard, Elevating, Mauve, Chrome (16")	1	1
29	B,E	1129606	Footboard, Elevating, Mauve, Grey (16")	1	1
29	B,D	1072303	Footboard, Elevating, Grotto, Chrome (16")	1	1
29	B,E	1129613	Footboard, Elevating, Grotto, Grey (16")	1	1
29	B,D	1072304	Footboard, Elevating, Blueberry, Chrome (16")	1	1
29	B,E	1129614	Footboard, Elevating, Blueberry, Grey (16")	1	1
29	B,D	1072297	Footboard, Elevating, Mauve, Chrome (18")	1	1
29	B,E	1129617	Footboard, Elevating, Mauve, Grey (18")	1	1
29	B,D	1072305	Footboard, Elevating, Grotto, Chrome (18")	1	1
29	B,E	1129615	Footboard, Elevating, Grotto, Grey (18")	1	1
29	B,D	1072306	Footboard, Elevating, Blueberry, Chrome (18")	1	1
29	B,E	1129616	Footboard, Elevating, Blueberry, Grey (18")	1	1

NOTE: A - Includes items 7, 14 and 15.  
 B - Includes items 1-28.  
 C - Item no longer available, refer to Part #1108498.  
 D - Before 1/24/05. For use on chairs with "Chrome Frame Finish."  
 E - After 1/23/05. For use on chairs with "Grey Frame Finish."  
 F - No Longer Available individually. See Category titled "Lateral Trunk Support" for Kit # HRA20K as replacement.