


# Pivot/Slide Tube w/ Footplate

THE PRODUCTION OF THE HTR IS CURRENTLY SUSPENDED AND SOME SERVICE PARTS MAY NOT BE AVAILABLE TO ORDER.



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		7500004-NLA	NLA - Pivot/Slide Tube w/ Plug Button	1	1
2		1027097	Plug Button, Black (3/4")	1	1
3		1006981	Spring, Tension	1	1
4		1015178	Footplate, Composite, Black (6" W x 7-53/64" L)	1	1
5		1016346	Cap, End, Tube, Black (7/8")	1	1
6		40120X001	Screw, Phillips Pan Head (#8-32 x 3/8")	1	2
7		1017224-NLA	NLA - Pivot/Slide Tube w/ Composite Footplate Assembly, Chrome - Right	1	1
7		1017225-NLA	NLA - Pivot/Slide Tube w/ Composite Footplate Assembly, Chrome - Left	1	1
	A	1027000	Heel Loop - Large - Right/Left	1	1
	A	1021421-NLA	NLA - Screw, Phillips Pan Head (1/4-20 x 2-1/4")	1	1
	A	1006477-NLA	NLA - Spacer (1/4 x 11/32 x 1-3/4")	1	1
	A	1025195	Locknut (1/4-20)	1	1

NOTE: A - Not Shown.