


Headrest Extension

THE PRODUCTION OF THE HTR IS CURRENTLY SUSPENDED AND SOME SERVICE PARTS MAY NOT BE AVAILABLE TO ORDER.



Headrest Extension

THE PRODUCTION OF THE HTR IS CURRENTLY SUSPENDED AND SOME SERVICE PARTS MAY NOT BE AVAILABLE TO ORDER.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		HRA18-NLA	NLA - Pad, Headrest	1	1
2		450106-25-NLA	Cover, Headrest Extension, Mauve	1	1
2		450106-26-NLA	Cover, Headrest Extension, Grotto	1	1
2		450106-27-NLA	Cover, Headrest Extension, Blueberry	1	1
3		1108143	Spacer, Nylon (3/16 x 5/16 x 5/16")	1	4
4		1013218	Washer (#10 Flat)	1	6
5	B	QC47A	Plate, Headwing - Chrome	1	1
5	C	1129545-NLA	Plate, Headwing - Grey	1	1
6	A	1108131-NLA	NLA - Screw, Phillips Pan Head (#10-24 x 1")	1	4
6		1108193-NLA	NLA - Screw, Phillips Pan Head (#10-24 x 3/4")	1	4
7		FB883	Knob, Adjustment	1	1
8	B	HRA1A	Bar, Slider, Chrome	1	1
8	C	1129514-NLA	Bar, Slider, Grey	1	1
9	D	1121008-NLA	NLA - Screw (#10-24 x 7/8")	1	2
10	B	46006925-NLA	Headrest Extension Assembly, Mauve, Chrome	1	1
10	C	1129620	Headrest Extension Assembly, Mauve, Grey	1	1
10	B	46006926-NLA	NLA - Headrest Extension Assembly, Grotto, Chrome	1	1
10	C	1129621	Headrest Extension Assembly, Grotto, Grey	1	1
10	B	46006927-NLA	Headrest Extension Assembly, Blueberry, Chrome	1	1
10	C	1129622	Headrest Extension Assembly, Blueberry, Grey	1	1

NOTE: A - Item no longer available, refer to Part #1108193.
 B - Before 1/24/05. For use on chairs with "Chrome Frame Finish."
 C - After 1/23/05. For use on chairs with "Grey Frame Finish."
 D - No Longer Available individually. See Category titled "Lateral Trunk Support" for Kit # HRA20K as replacement.