


# Mattress Assembly - APM

For Units Manufactured After Serial Number 99G...



## Mattress Assembly - APM

### For Units Manufactured After Serial Number 99G...

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		9830002-NLA	NLA - Shell, Mattress	1	1
2		9830003-NLA	NLA - Cushion, Air	1	20
3		9830004-NLA	NLA - Manifold, APM	1	1
4		9650062-NLA	NLA - Pad, Therapy	1	1
5		9830007-NLA	NLA - Hose Set, Mattress	1	1
6	B	BB9652050-NLA	NLA - Mattress Assembly	1	1
	A	9570020-NLA	NLA - Mat, Transport	1	1
	A	BB9830006-NLA	NLA - Carton	1	1
	A	9170137	NLA - Bag, Duffel	1	1
NOTE: A - Not Shown B - Includes items 1-5, Mat, Carton, and Duffle Bag.					