


Mattress Assembly - LTM



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		9720025-NLA	NLA - Pad, Therapy	1	1
2		9720027-NLA	NLA - Mattress Cover	1	1
3		9720026-NLA	NLA - Tube, Air	1	8
4		9720024	Mattress Lower, Turn-Q	1	1
5		833	Elbow, Nylon (5/8")	1	4
6		268	Hose, Spiral (100' Roll)	1	1
7		BB9722039-NLA	NLA - Tee, Side Ventilated (5/8")	1	1
7		9722039-NLA	NLA - Tee, Side Ventilated (5/8")	1	1
8		1101069-NLA	NLA - Tee (5/8")	1	5
9		9652039	Plate, Coupler, CPR	1	1
10		9722036-NLA	NLA - Hose Set Assembly (Includes: Right, Left, and Center Hoses)	1	1
	A	BB9722002-NLA	NLA - LTM Mattress Assembly, Complete	1	1

NOTE: A - Not Shown