

**APM
CG™ 8000
LTM
MNS Series
Mattress**



Yes, you can.

Parts Catalog Usage Guide

The information contained in this document is subject to change without notice. This catalog contains service parts for Form No. 96-143. If you have any questions regarding this catalog, please contact one of our Customer Service Representatives.

For Customer Service Support Call:

1-800-333-6900

For Faxed In Orders Call:

1-800-678-4682 or

1-440-366-0709

For Technical Service Call:

1-800-832-4707

Literature Can Be Found at:

www.invacare.com or

1-800-828-6272

SYMBOL DEFINITION

Throughout this catalog you will notice letters to draw your attention to specific footnotes which are located at the bottom of the page or below specified charts.

Examples:

A,B,C etc. = Special Comments

NOTE: = General information pertaining to that specific page

NLA = No Longer Available

HARDWARE MEASUREMENTS

The measurements for the hardware will be listed as follows:

Screws, Bolts Diameter X Thread Count/Pitch X Length

Washers Inside Diameter X Outside Diameter X Thickness

Spacers Inside Diameter X Outside Diameter X Length

For your convenience, we have supplied an area for you to list your Invacare Representatives.

Name	Phone Numbers


©2019 Invacare Corporation, Elyria, Ohio. Invacare, Yes, you can. and the Medallion Design are registered trademarks; CareGuard is a Trademark of the Invacare Corporation.

TABLE of CONTENTS

Description	Page Number
Careguard 8000, Mattress, and Pump	10
Mattress Assembly - APM (After Serial No. 99G)	7
Mattress Assembly - APM (Before Serial No. 99H)	6
Mattress Assembly - LTM	9
Miscellaneous Accessories - APM and LTM	12
MNS400 APM System Service Parts	13
MNS500 Low Air Loss System Service Parts	15
MNS600 System Service Parts	17
Packaging and Literature - Mattress	3
Packaging and Literature - Power Unit	4
Power Unit Assembly	5


Packaging and Literature - Mattress

MATTRESS



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		BB9722005-NLA	NLA - Duffel Bag Assembly - LTM	1	1
1		9652009-NLA	NLA - Duffel Bag Assembly - APM	1	1
2		BB9170153-NLA	NLA - Liner, Carton	1	1
3		9170098-NLA	NLA - Carton, Mattress	1	1
NOTE: When adjusting or replacing parts, please refer to your Owner's Manual for assistance.					


Packaging and Literature - Power Unit



NOTE: When adjusting or replacing parts, please refer to your Owner's Manual for assistance.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		9570023-NLA	NLA - Carton	1	1
2		1101068-NLA	NLA - Bag, Poly	1	1
3		9652026-NLA	NLA - Cap, End, Foam Molded, APM - Right	1	1
3		9652042-NLA	NLA - Cap, End, Foam Molded, LTM - Right	1	1
4		9652027-NLA	NLA - End, Cap, Foam Molded, APM - Left	1	1
4		9652043-NLA	NLA - Cap, End, Foam Molded, LTM - Left	1	1
5		9650040	Owner's Manual - APM	1	1
5		9720008-NLA	NLA - Owner's Manual - LTM	1	1

Power Unit Assembly




Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		BB9650033 -NLA	NLA - Cover, Front	1	1
2		817	O-Ring (5mm x 2mm)	1	4
3		9540110-NLA	NLA - Plate, Hook	1	2
4		9540111-NLA	NLA - Hook, Narrow	1	2
5		9540112-NLA	NLA - Hook, Wide	1	2
6		9650029-NLA	NLA - Gasket, Hook Plate	1	2
7		9650019	Element, Air Filter	1	1
8	A	845-NLA	NLA - V-Ring, CPR	1	3
9	A	424-NLA	NLA - O-Ring (15mm ID x 3mm)	1	3
	B,C	9652001-NLA	NLA - Power Unit, APM Boxed	1	1
	B,C	9653001-NLA	NLA - Chassis, APM Boxed	1	1
	B,C	BB9722001-NLA	NLA - Power Unit, LTM Boxed	1	1

NOTE: A - For LTM units only and for APM units manufactured before serial number 99H.
 B - Not Shown
 C - For a complete boxed APM Power Unit; before Serial Number 99H.. order part no. 9652001, after Serial Number 99G.. is not available.

Mattress Assembly - APM

For Units Manufactured Before Serial Number 99H...




Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		9650046-NLA	NLA - Core, Mattress	1	4
2		9650009-NLA	NLA - Mattress Case	1	1
3		9650030	Hook and Loop (2" W x 5-1/2" H)	1	4
4		9650031-NLA	NLA - Hook and Loop (2 x 2")	1	4
5		9652011-NLA	NLA - Hose, Channel A	1	1
6		9652012-NLA	NLA - Hose, Channel B	1	1
7		9652013-NLA	NLA - Hose, Channel C	1	1
8		9652039	Plate, Coupler, CPR	1	1
9		BB9650004-NLA	NLA - Pad, Therapy	1	1
	A	9652002-NLA	NLA - Boxed Mattress Assembly	1	1

NOTE: A - Not Shown.

Mattress Assembly - APM

For Units Manufactured After Serial Number 99G...




Mattress Assembly - APM

For Units Manufactured After Serial Number 99G...

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		9830002-NLA	NLA - Shell, Mattress	1	1
2		9830003-NLA	NLA - Cushion, Air	1	20
3		9830004-NLA	NLA - Manifold, APM	1	1
4		9650062-NLA	NLA - Pad, Therapy	1	1
5		9830007-NLA	NLA - Hose Set, Mattress	1	1
6	B	BB9652050-NLA	NLA - Mattress Assembly	1	1
	A	9570020-NLA	NLA - Mat, Transport	1	1
	A	BB9830006-NLA	NLA - Carton	1	1
	A	9170137	NLA - Bag, Duffel	1	1
NOTE: A - Not Shown B - Includes items 1-5, Mat, Carton, and Duffel Bag.					


Mattress Assembly - LTM



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		9720025-NLA	NLA - Pad, Therapy	1	1
2		9720027-NLA	NLA - Mattress Cover	1	1
3		9720026-NLA	NLA - Tube, Air	1	8
4		9720024	Mattress Lower, Turn-Q	1	1
5		833	Elbow, Nylon (5/8")	1	4
6		268	Hose, Spiral (100' Roll)	1	1
7		BB9722039-NLA	NLA - Tee, Side Ventilated (5/8")	1	1
7		9722039-NLA	NLA - Tee, Side Ventilated (5/8")	1	1
8		1101069-NLA	NLA - Tee (5/8")	1	5
9		9652039	Plate, Coupler, CPR	1	1
10		9722036-NLA	NLA - Hose Set Assembly (Includes: Right, Left, and Center Hoses)	1	1
	A	BB9722002-NLA	NLA - LTM Mattress Assembly, Complete	1	1

NOTE: A - Not Shown

Careguard 8000, Mattress, and Pump



Careguard 8000, Mattress, and Pump

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		9830002-NLA	NLA - Shell, Mattress	1	1
2		9830003-NLA	NLA - Cushion, Air	1	20
3		9830004-NLA	NLA - Manifold, APM	1	1
4		9830005-NLA	NLA - Therapy Pad	1	1
5		9830007-NLA	NLA - Hose Set, Mattress	1	1
5		9830008-NLA	NLA - Hose Set - Pump	1	1
6		877-NLA	NLA - Pump	1	1
7	B	BB9652050-NLA	NLA - Mattress Assembly	1	1
	A	9570020-NLA	NLA - Mat, Transport	1	1
	A	9650040	Owner's Manual - APM	1	1
	A	BB9830006-NLA	NLA - Carton	1	1
	A	9170137	NLA - Bag, Duffel	1	1
NOTE: A - Not Shown B - Includes items 1-5, Mat, Manual, Duffel Bag					

Miscellaneous Accessories - APM and LTM

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
		9652031-NLA	NLA - Mattress Core Set	4	1
		9652032-NLA	NLA - Hose Set	2	1
		BB9612031-NLA	NLA - Extension Set, Side Rail - Full	2	1
		BB9612030-NLA	NLA - Extension Set, Side Rail - Half	2	1
		9540048	NLA - Sheet, Draw (6 x 3")	1	1
		9570020-NLA	NLA - Mat, Transport	1	1
		BB730-NLA	NLA - Pads, Incontinence	72	1

MNS400 APM System

Mattress Covers are not shown.



MNS400 APM System


Mattress Covers are not shown.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1135275-NLA	NLA - Cover, Mattress Top MNS400 - Narrow (STZ)	1	1
1	A	1135276-NLA	NLA - Cover, Mattress Top MNS400 - Standard (STZ)	1	1
1	A	1135277-NLA	NLA - Cover, Mattress Top MNS400 - Bariatric (STZ)	1	1
2	A	1132039-NLA	NLA - Cover, Mattress Top MNS400 - Narrow (LTZ)	1	1
2	A	1132040-NLA	NLA - Cover, Mattress Top MNS400 - Standard (LTZ)	1	1
2	A	1132041-NLA	NLA - Cover, Mattress Top MNS400 - Bariatric (LTZ)	1	1
3		1141540-NLA	NLA - Connector, Male CPC (1/4" NPT)	1	1
4		1141539-NLA	NLA - Connector, Female CPC 1/4" NPT	1	1
5		1141543-NLA	NLA - Fuse (750 M AMP)	5	1
6		1141587-NLA	NLA - Package, Vibration Mount	4	1
7		1141542-NLA	NLA - Muffler, Brass (1/8" NPT)	1	1
8		1141533-NLA	NLA - Package, Solenoid Valve	2	1
9		1141541-NLA	NLA - Fitting, Elbow (1/4" Barb x 1/8" NPT)	2	1
10		1141536-NLA	NLA - Compressor, Apollo	1	1
11		1143685-NLA	NLA - Tubing, Braided Inner (25' L)	1	1
12		1143684-NLA	NLA - Tubing, Vinyl Double (25' L)	1	1
13		1146110-NLA	NLA - Bladder, Air Sack, MNS400 Narrow (32" Long)	1	1
13		1143680-NLA	NLA - Bladder, Air Sack, MNS400 & MNS400-E Standard (35" Long)	1	1
13		1146111-NLA	NLA - Bladder, Air Sack, MNS400 Bariatric (42" Long)	1	1

NOTE: A - Not shown

MNS500 Low Air Loss System

Mattress Covers are not shown.




MNS500 Low Air Loss System

Mattress Covers are not shown.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1135278-NLA	NLA - Cover, Mattress Top MNS500 - Narrow (STZ)	1	1
1	A	1135279-NLA	NLA - Cover, Mattress Top MNS500 - Standard (STZ)	1	1
1	A	1135280-NLA	NLA - Cover, Mattress Top MNS500 - Bariatric (STZ)	1	1
2	A	1132042-NLA	NLA - Cover, Mattress Top MNS500 - Narrow (LTZ)	1	1
2	A	1132043-NLA	NLA - Cover, Mattress Top MNS500 - Standard (LTZ)	1	1
2	A	1132044-NLA	NLA - Cover, Mattress Top MNS500 - Bariatric (LTZ)	1	1
3		1141540-NLA	NLA - Connector, Male CPC (1/4" NPT)	1	1
4		1141539-NLA	NLA - Connector, Female CPC 1/4" NPT	1	1
5		1141543-NLA	NLA - Fuse (750 M AMP)	5	1
6		1141587-NLA	NLA - Package, Vibration Mount	4	1
7	B	1141522-NLA	NLA - Package, Thomas Filter	10	1
8		1141542-NLA	NLA - Muffler, Brass (1/8" NPT)	1	1
9		1141533-NLA	NLA - Package, Solenoid Valve	2	1
10		1141561-NLA	NLA - Insert, Coupling (1/2" NPT)	1	1
11	C	See Note	See Footnote	1	1
12		1146109-NLA	NLA - Bladder, Air Sack, MNS5000 Narrow (32" Long)	1	1
12		1143681-NLA	NLA - Bladder, Air Sack, MNS5000 Standard (35" Long)	1	1
12		1143682-NLA	NLA - Bladder, Air Sack, MNS5000 Bariatric (42" Long)	1	1
NOTE: A - Not shown B - Photo shows 1(one) filter, part # contains a quantity of 10 (ten) Filters. C - No longer available as a service part in the field.					

MNS600 System

Mattress Covers are not shown.



MNS600 System

Mattress Covers are not shown.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1135281-NLA	NLA - Cover, Mattress Top MNS600 - Narrow (STZ)	1	1
1	A	1135282-NLA	NLA - Cover, Mattress Top MNS600 - Standard (STZ)	1	1
1	A	1135283-NLA	NLA - Cover, Mattress Top MNS600 - Bariatric (STZ)	1	1
2	A	1132045-NLA	NLA - Cover, Mattress Top MNS600 - Narrow (LTZ)	1	1
2	A	1132046-NLA	NLA - Cover, Mattress Top MNS600 - Standard (LTZ)	1	1
2	A	1132047-NLA	NLA - Cover, Mattress Top MNS600 - Bariatric (LTZ)	1	1
3		1141540-NLA	NLA - Connector, Male CPC (1/4" NPT)	1	1
4		1141539-NLA	NLA - Connector, Female CPC 1/4" NPT	1	1
5		1141543-NLA	NLA - Fuse (750 M AMP)	5	1
6		1141587-NLA	NLA - Package, Vibration Mount	4	1
7		1141542-NLA	NLA - Muffler, Brass (1/8" NPT)	1	1
8		1141533-NLA	NLA - Package, Solenoid Valve	2	1
9		1141561-NLA	NLA - Insert, Coupling (1/2" NPT)	1	1
10		1141541-NLA	NLA - Fitting, Elbow (1/4" Barb x 1/8" NPT)	2	1
11		1141572-NLA	NLA - Solenoid Manifold	1	1
12		1141574-NLA	NLA - Compressor - Thomas	1	1
13		1143684-NLA	NLA - Tubing, Vinyl Double (25' L)	1	1
14		1146106-NLA	NLA - Bladder, Air Sack, MNS600 Narrow (32" Long)	1	1
14		1146107-NLA	NLA - Bladder, Air Sack, MNS600 Standard (35" Long)	1	1
14		1146108-NLA	NLA - Bladder, Air Sack, MNS600 Bariatric (42" Long)	1	1
NOTE: A - Not shown					



Invacare Corporation

www.invacare.com

One Invacare Way
Elyria, Oh USA
44036-2125
800-333-6900

Technical Support
800-832-4707
Florida Technical Support
407-321-5630

Customer Service
800-333-6900

©2019 Invacare Corporation, Elyria, Ohio.
All rights reserved. Trademarks are
identified by the symbols ™ and ®. All
trademarks are owned by or licensed to
Invacare Corporation unless otherwise
noted. (Revised 04/10/19)

Form No. 96-143