

# Footrest & Legrest Accessories

## Footrest and Legrest Accessories

Part Number	Description	Quantity Per Package	Quantity Per Assembly
1036099-NLA See Note(s) A	Heel Loops w/o Hardware (Pair)	2	1
1036100 See Note(s) B	Heel Loops w/o Hardware (Pair)	2	1
1035905 See Note(s) C	Heel Loop w/o Hardware (Pair) - Small	2	1
1036096 See Note(s) D	Hardware, Heel Loop (Pair)	2	1
1035926 See Note(s) E	Hardware, Heel Loops Aluminum Footplate (Pair)	2	1
1143050	Non-Removable Device	1	1
NOTE: A - Use on a 10" Wide Chair. B - J932A, J932A-14, J933A and J934A Front Riggings. C - J934A-14 Front Riggings. D - Use with Heel Loop #1036099 and 1036100. E - Use with Heel Loop #105905.			