

Invacare®  
Pronto™ M51™  
Power  
Wheelchairs



Invacare® Pronto™ M51™ Power Wheelchairs



Yes, you can.™



Yes, you

Presented in a stylish, compact size, the Invacare® Pronto™ M51™ power wheelchair with SureStep™ has true center-wheel drive performance for exceptional maneuverability and intuitive driving. Featuring a medium back van seat and innovative SureStep technology, the Pronto M51 provides a smooth, less jarring ride over transitions and thresholds up to two inches in height while maintaining stability. The Pronto M51 is attractively designed for optimum maneuverability and accessibility.



## SURE STEP

### SureStep™ Technology

SureStep technology provides optimum traction and stability as the Pronto M51 drives over transitions and uneven ground up to two-inches in height. The powerful motors lift the front casters to assist with changes in height for a smooth, less jarring ride over rough surfaces. SureStep technology allows six wheels to remain in contact with the ground to maximize performance, improve stability and enhance comfort for a more controllable power chair.



Front casters lift up to accommodate transitions up to two-inches for a smooth and controllable ride.

SureStep design flexes to maintain drive-wheel traction and pulls the chair over obstacles and transitions.



Six wheels remain on the ground for maximum stability, maneuverability and accessibility.

### FEATURES

- SureStep technology allows for smooth driving over thresholds and transitions
- Stylish design with choice of Candy Red or Deep Blue shroud colors
- Compact size for the tightest overall turning radius possible
- True center-wheel-drive system for intuitive driving
- Six wheels on the ground provide optimum stability
- Fully programmable MKIV™-R11 electronics for maximum function and precise driver control
- Adjustable joystick mount
- Medium Back van seats and with flip-up adjustable height and width arms
- Powerful 350-watt, high-torque motors for reliable performance
- Regenerative braking and automatic motor locks for dependable stopping
- 3-amp on-board charger for convenience
- 300 lb. weight capacity

### Specifications

Overall dimensions	
Lengths (with footboard folded)	34"
Lengths (with footboard down)	38" - 40"
Width (excluding joystick)	24"
Heights	35" - 37"
Seat dimensions	
Width	18"
Depth	16"-18"
Weight	
With batteries	205 lb.
Without batteries	155 lb.
Wheel and caster sizes	
Front casters	6" x 1.3"
Drive wheels (flat-free)	10.5" x 3.5"
Rear casters	6" x 1.3"
Maximum speed	3 mph
Incline capability	9°
Weight capacity	300 lb.
Turning radius	
Front (with footboard down)	19.5"
Rear	19.5"
Ground clearance	3.8"
Warranties	
Base frame	60 months
Seat frame	60 months
Electronics	12 months
Motor/gearbox	12 months

Invacare Corp  
www.invacare.c

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

Invacare is a registered trademark and MKIV, Pronto, SureStep you can. are trademarks of Invacare Corporation. Specifications are subject to change without notice.

© 2002

Invacare Corporation  
Form No. 02-066  
Rev. 10/02



"Invacare says 'Yes' and I couldn't have better myself."

*Arnold Hatanaka*