



# WalkLite Walkers

Yes, you can.®



NEW

The Invacare® WalkLite™ Walker revolutionized the world of walkers by addressing the needs of providers, clinicians and consumers at an affordable price. The new WalkLite line of walkers covers a wide range of consumer heights from Junior to Tall Adult. The new WalkLites will continue to feature all the benefits of the original WalkLite – a comfortable flip-up seat, durable front 5" fixed wheels, rear spring-loaded SureGlide™ brakes, a lightweight aluminum frame with flexible backrest and ergonomic dual-paddle folding mechanisms.



**Yes, you can.®**

**Invacare Corporation**  
www.invacare.com

**USA**

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**

570 Matheson Blvd. E.  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324



Model no. 6300-JRA



Model no. 6300-ATA



*Folds into a compact unit for storage. Built-in seat handle allows for easy transport.*

**Specifications**

	6300-JRA	6300-ATA
Depth opened at base	23"	23"
Depth folded	7"	7"
Width inside handgrips	18"	18"
Width inside base of legs	20"	24.5"
Width overall	27"	27"
Height adjustment min.-max.	29.25" - 36.50"	32" - 39"
Approx. Client height* min.-max.	4'11" - 6'0"	5'5" - 6'5"
Seat-to-floor height min.-max.	16.625" - 23.625"	19" - 25.75"
Seat width x depth	17" x 9.50"	17" x 9.50"
Wheels	5" fixed	5" fixed
Product weight	11 lb.	11 lb.
Weight capacity	300 lb.	300 lb.
Quantity	1 per carton	1 per carton
Limited warranty	lifetime frame	lifetime frame
other components	1 year	1 year

\*Client height is approximation. Each individual should always consult with his/her physician or therapist to determine proper adjustment.

**Features**

- Lightweight aluminum frame
- Flip-up seat with built-in handle
- Ergonomic dual-paddle folding mechanisms
- Folds into a compact unit for storage and transport
- Flexible backrest for comfort
- Rear spring-loaded SureGlide brakes
- Adjustable seat-to-floor height



**Basket and Tray Accessories**

- Basket (model no. 65110 for 6300-ATA, 65109 not shown for 6300-JRA)
  - 10 lb. weight capacity
  - Position basket underneath seat to secure valuables
- Tray (model no. 65408)
  - For food and beverage placement
  - Tray fits on both 6300-JRA and 6300-ATA



*Flip-up seat allows user to step inside the frame.*



*Rear spring-loaded SureGlide brakes provide additional security when seated.\*\**

\*\* The rear SureGlide brake's engagement may vary depending upon the weight distributed on the seat and the positioning of the user on the seat.

Invacare, the Blue Medallion Design and Yes, you can., are registered trademarks and WalkLite and SureGlide are trademarks of Invacare Corporation. Specifications are subject to change without notification.

Patents: 5,529,425 D480,666 D480,995 and D483,702. Patents pending.

© 2005

Invacare Corporation

Form No. 03-102

Rev. 10/05



*"Invacare says 'Yes, you can,' and I couldn't have said it better myself."*